

auto-Lancet™

“Custom”



Enduring Quality with a Modern Design

The incomparable **auto-Lancet™** is the lancing device that has set the standard for quality, comfort, and adjustability since 1983. This re-usable lancing device helps patients obtain capillary blood samples for testing. The auto-Lancet “Custom” adds a new look and size to the famous line of re-usable lancing devices from Palco Labs. Meter and kit manufacturers have a wide selection of barrel colors to match any design, in addition to a customized, debossed logo. The new design is even more compact and comfortable than the traditional auto-Lancet “Classic.”

Quality... *that set the standard for re-usable lancing devices and became the choice for glucose meter manufacturers worldwide*

- Same durable, high quality materials
- Compatible with standard industry lancets
- Same precision engineering
- Same Lifetime Warranty

Advanced Design... *in the tradition of Palco engineering excellence*

- More comfortable to grip and activate
- Soft, elastomer barrel in wide selection of colors for customized logo
- Ergonomic adjustable tip
- More compact size in a pen shape
- Faster action and **even more** linear strike than standard-sized auto-Lancet “Classic”

Enhanced Performance... *meeting the needs of today's diabetic market*

- A more linear impact for minimal cell damage and consistent blood sample
- A faster impact to further reduce the sensation and expectation of pain
- Additional depth settings to suit individualized needs for comfort and blood sample

PALCOLABS

Specifications

auto-Lancet "Custom" Specifications:

Length: 4.905 in. (12.45 cm)

Barrel Width: 0.475 in. (1.206 cm)

Materials: Polycarbonate (Threaded Center), Acrylonitrile Butadiene Styrene (Hub, Sleeve), Acetal (Carrier)

Lancet Compatibility:

Industry standard, including Palco Labs EZ-Lets, Ultra-Fine* II, Monolet*, Unilet* GP, Microlet*, OneTouch*, Ultra TLC*, SoftTouch* [Roche Australia], SoftTouch* [Roche Accu-Chek], Gentle-Let*, CleanLet*, Omnican*, Qualitilet*, and Freestyle*.

HCPC Code: A4259

Warranty:

The auto-Lancet is guaranteed against manufacturing defects for the lifetime of the product.

Ordering Information

Palco Order Hotline: 1-800-346-4488

Order Fax: 831-430-1610

E-mail: orders@palcolabs.com

Please note that not all products are available in all markets. For full details, please contact your local reseller or Palco Labs.

Special Orders

For special configurations and high volume orders, contact Palco Labs or send e-mail to config@palcolabs.com.

auto-Lancet is a trademark of Palco Labs, Inc.

U.S. Patent No. 5,613,978;

Germany Patent No. 19,781,810.2;

Hungary Patent No. 220,178;

Japan Patent No. 2,815,706.

All rights reserved.

*Other trademarks are property of their respective holders.

PALCOLABS

360 El Pueblo Road, Suite 102

Scotts Valley, CA • 95066 USA

Ph: 1-831-430-1600 • Toll-free: 1-800-346-4488

e-mail: info@palcolabs.com

web: www.palcolabs.com

Specifications are subject to change without notice.

© 2006 Palco Labs, Inc. All rights reserved.

1001-10302-100