

# EZ-Lance<sup>®</sup>

## The Ultimate Safety Lancet

The revolutionary new design of the **EZ-Lance<sup>®</sup>** maximizes comfort, safety, and ease of use. Every feature has been designed with users in mind, from the extra-large tab for easy removal to the ergonomic thumb-groove for the utmost in comfort. The high-speed, linear-tracking lance is highly advanced technology designed to help minimize pain. The internal breakaway bar is designed to prevent reuse, while the safety indicator window on the bottom ensures that you know whether the **EZ-Lance** is used or unused.

Don't settle for less than the best - demand the **EZ-Lance!**

From the people at Palco Labs, makers of world-class medical device products.

Ideal for use in  
hospitals and labs



### Key Features

#### **Safety**

- Single-use to eliminate cross-contamination
- Not re-cockable - maximum safety feature prevents reuse
- Safety indicator shows clearly if device has been used
- No needle exposure for maximum safety; needle is covered at all times before and after lance

#### **Comfort**

- Ergonomic design - easy to hold; large tab is easy to remove
- High-speed lance - helps minimize pain
- Linear tracking - minimizes needle vibration for a cleaner lance
- Special bounce-back restraint - single lance only, no chance for a second lance

#### **Ease of Use**

- Advanced design - minimal user instructions
- Pre-cocked - simplifies use; no extra step needed

### Key Benefits

- Four different versions available
- Maximum safety - impossible to re-cock, lance is always concealed except at the moment of lancing
- Ultimate ease of use - remove tab, apply to skin, and press button to fire
- Ergonomic design for maximum comfort
- Advanced lancing technology to minimize pain

## Specifications & Packaging Configurations

### Ordering Information

Palco Order Hotline: 1-800-346-4488

Order Fax: 831-430-1610

E-mail: [info@palcolabs.com](mailto:info@palcolabs.com)

Please note that not all products are available in all markets. For full details, please contact your local reseller or Palco Labs.

### Special Orders

For special configurations and high volume orders, contact Palco Labs or send e-mail to [info@palcolabs.com](mailto:info@palcolabs.com).

#### Super Soft Model

Color: White (body) and yellow (tab)

Needle size: 26 gauge

Point configuration: Tri-bevel

Depth of Penetration: 1.8 mm



**PN/REF: ELY100** (100 ct. box)

UPC Code: 0-25294-04814-9

#### Soft Model

Color: White (body) and green (tab)

Needle size: 21 gauge

Point configuration: Tri-bevel

Depth of Penetration: 1.8 mm



**PN/REF: ELG100** (100 ct. box)

UPC Code: 0-25294-04808-8

#### Normal Model

Color: White (body) and orange (tab)

Needle size: 21 gauge

Point configuration: Tri-bevel

Depth of Penetration: 2.2mm



**PN/REF: ELO100** (100 ct. box)

UPC Code: 0-25294-04810-1

#### Ultra Model

Color: White (body) and red (tab)

Needle size: 21 gauge

Point configuration: Microblade

Depth of Penetration: 2.2mm



**PN/REF: ELR100** (100 ct. box)

UPC Code: 0-25294-04812-5

#### All Models

Length (w/tab): 1.87 in. (47.49 mm)

Largest width: 0.79 in. (20.06 mm)

Height: 0.54 in. (13.71 mm)

Materials: Polyethylene/polypropylene/copolymer (needle carrier), stainless steel (needle), polystyrene (case) - **Latex free**

**HCPC Code (all models):** A4259

#### Warranty (all models):

Palco Labs product carries a one (1) year warranty against manufacturer's defects in materials and workmanship.

EZ-Lance is a registered trademark of Palco Labs, Inc.  
U.S. and international patents pending.

Registered Community Design No. 0001139020001

Japanese Design Registration No. 1,214,580

All Rights Reserved.

## PALCOLABS

360 El Pueblo Road, Suite 102  
Scotts Valley CA • 95066, USA

Ph: 1-831-430-1600 • Toll-free: 1-800-346-4488

e-mail: [info@palcolabs.com](mailto:info@palcolabs.com)

web: [www.palcolabs.com](http://www.palcolabs.com)

Specifications listed are preliminary. Final specifications subject to change without notice.

© 2006, Palco Labs, Inc.

1001-60001-006