


**PALCOLABS**

# EZ Lancets™

---



## Premium Blood Lancets

*Experience the best in blood lancets from Palco Labs, a global leader in the design and manufacture of medical device products focused on reducing pain. EZ-Lancets are precision engineered for the finest in needle edges, reducing discomfort and encouraging the healing process. EZ-Lancets fit most re-usable finger lancing devices, and are available in three sizes: 21 gauge, 26 gauge and 30 gauge.*

*Rely on Palco Labs to deliver the highest quality blood lancets for capillary sampling requirements!*

---

### Key Features

#### **Safety**

- Sterile blood lancets
- Single use to eliminate cross-contamination
- Twist cap to expose needle
- Cap may be used to receive used needle, to avoid blood contamination

#### **Precision Engineered**

- Finely-ground needle edges for cleaner lance
- Fine points for least amount of capillary damage
- Consistent quality

### Key Benefits

- Three sizes available: 21 gauge, 26 gauge and 30 gauge
- Available in 100 count and 200 count boxes
- Ultimate ease of use - insert in device and remove needle cap to use
- Highest quality for maximum comfort and least amount of cell damage

## Ordering Information

Palco Order Hotline: 1-800-346-4488  
Order Fax: 831-430-1610  
E-mail: [info@palcolabs.com](mailto:info@palcolabs.com)

Please note that not all products are available in all markets. For full details, please contact your local reseller or Palco Labs.

## Special Orders

For special configurations and high volume orders, contact Palco Labs or send e-mail to [info@palcolabs.com](mailto:info@palcolabs.com).

## Specifications & Packaging Configurations

### 30 Gauge

Needle size: 30 gauge  
Point configuration: Tri-bevel

**30 gauge**

**PN/REF: ETP100** (100 ct. box)

UPC Code: TBD

**PN/REF: ETP200** (200 ct. box)

UPC Code: TBD

### 26 Gauge

Needle size: 26 gauge  
Point configuration: Tri-bevel

**26 gauge**

**PN/REF: ETY100** (100 ct. box)

UPC Code: TBD

**PN/REF: ETY200** (200 ct. box)

UPC Code: TBD

### 21 Gauge

Needle size: 21 gauge  
Point configuration: Tri-bevel

**21 gauge**

**PN/REF: ETO100** (100 ct. box)

UPC Code: TBD

**PN/REF: ETO200** (200 ct. box)

UPC Code: TBD

### All Models

Length (w/tab): 1.5 in. (38.1 mm)

Diameter: 0.25 in. (6.35 mm)

Materials: Polypropylene (needle carrier), stainless steel (needle), Latex free

### For use with most industry-standard finger-lancing devices, including:

auto-Lancet™, SoftTouch™, Qwik-Let™, Penlet® Plus, Penlet® II,  
MicroLet™ Vaculance, Dialet™, OneTouch UltraSoft™, Monojector™, Freestyle®,  
B-D • Medisense®, Microlet™, Autolet®

### Warranty (all models):

This Palco Labs product carries a one (1) year warranty against manufacturer's defects in materials and workmanship.

EZ-Lancets is a trademark of Palco Labs, Inc.  
U.S. and international patents pending.  
All Rights Reserved.

## PALCOLABS

360 El Pueblo Road, Suite 102  
Scotts Valley CA • 95066, USA  
Ph: 1-831-430-1600 • Toll-free: 1-800-346-4488  
e-mail: [info@palcolabs.com](mailto:info@palcolabs.com)  
web: [www.palcolabs.com](http://www.palcolabs.com)

Specifications listed are preliminary. Final specifications subject to change without notice.  
© 2006, Palco Labs, Inc.