

auto-Lancet™

Enduring Quality

The Original Adjustable Tip

The incomparable **auto-Lancet™** is the lancing device that has set the standard for quality, comfort, and adjustability since 1983. This re-usable lancing device assists in obtaining capillary blood samples for testing. The **auto-Lancet** is constructed in the USA of durable plastics for precision lancing day after day, year after year. It's the only lancing device with a lifetime warranty.

The **auto-Lancet** - quality you can count on.
From the people at Palco Labs,
makers of world-class medical device
products and leaders in
lancing technology.



Key Features/Benefits

- Adjustable tip - Allows precise customization of lancet penetration
- High quality construction - Under lifetime warranty
- Made of durable plastics for years of use
- Convenient pen shape available in standard or mini size
- Made in the USA

Key Selling Points

- Available in standard or mini size
- Works with industry standard lancets
- #1 lancing device chosen by glucose meter manufacturers worldwide

PALCOLABS

Ordering Information

Palco Order Hotline: 1-800-346-4488
Order Fax: 831-430-1610
E-mail: orders@palcolabs.com

Please note that not all products are available in all markets. For full details, please contact your local reseller or Palco Labs.

Special Orders

For special configurations and high volume orders, contact Palco Labs or send e-mail to config@palcolabs.com.

auto-Lancet is a trademark of Palco Labs, Inc.
U.S. Patent No. 5,613,978;
Germany Patent No. 19,781,810.2;
Hungary Patent No. 220,178;
Japan Patent No. 2,815,706.
All rights reserved.

*Other trademarks are property of their respective holders.

PALCOLABS

360 El Pueblo Road, Suite 102
Scotts Valley, CA • 95066 USA
Ph: 1-831-430-1600 • Toll-free: 1-800-346-4488
e-mail: info@palcolabs.com
web: www.palcolabs.com

Specifications are subject to change without notice.
© 2005 Palco Labs, Inc. All rights reserved.

Specifications



auto-Lancet Standard (actual size)

Standard auto-Lancet Specifications:

Length: 4.905 in. (12.45 cm)
Barrel Width: 0.475 in. (1.206 cm)
Materials: Polycarbonate (Threaded Center), Acrylonitrile Butadiene Styrene (Hub, Sleeve), Acetal (Carrier)



auto-Lancet Mini (actual size)

Mini auto-Lancet Specifications:

Length: 3.595 in. (9.13 cm)
Barrel Width: 0.470 in. (1.19 cm)
Materials: Polycarbonate (Threaded Center), Acrylonitrile Butadiene Styrene (Hub, Sleeve), Acetal (Carrier)

Lancet Compatibility:

Industry standard, including Palco Labs EZ-Lets, Ultra-Fine* II, Monolet*, Unilet* GP, Microlet*, OneTouch*, Ultra TLC*, SoftTouch* [Roche Australia], SoftTouch* [Roche Accu-Chek], Gentle-Let*, CleanLet*, Omnican*, Qualitilet*, and Freestyle*.

HCPC Code: A4259

Warranty:

Palco Labs warrants that each new auto-Lancet shall be free of all manufacturing defects and guarantees normal performance for the lifetime of the owner.

Available Packaging Configurations

PN/REF: AL100K6AS

Contents: Clamshell package with one standard size auto-Lancet, tray filler card, instructions for use, and two 23G EZ-Lets
UPC Code: 0-25294-83100-0

PN/REF: AL100MK6AS

Contents: Clamshell package with one mini size auto-Lancet, tray filler card, instructions for use, and two 23G EZ-Lets
UPC Code: 0-25294-00100-7

PN/REF: AL100K6ASP

Contents: Soft Pack (poly bag) with one standard size auto-Lancet, tray filler card, instructions for use, and two 23G EZ-Lets
UPC Code: 0-25294-00601-9

PN/REF: AL100MK6ASP

Contents: Soft Pack (poly bag) with one mini size auto-Lancet, tray filler card, instructions for use, and two 23G EZ-Lets
UPC Code: 0-25294-00600-2

PN/REF: AL100K6AESP

Contents: Soft Pack (poly bag) with one standard size auto-Lancet, tray filler card, instructions for use, and NO EZ-Lets
UPC Code: 0-25294-00602-6

PN/REF: AL100K6AB

Contents: Soft Pack (poly bag) with 100 standard size auto-Lancets, 10 poly bags per carton.

PN/REF: AL100MK6AB

Contents: Soft Pack (poly bag) with 100 mini size auto-Lancets, 10 poly bags per carton.